



## TG-142 MADE EASY

Expedite superior QA, using an intuitive user-interface, specialized phantoms and centralized database to simplify over 30 TG-142 procedures



# PIPSPRO SOFTWARE

### ● COMPREHENSIVE IMAGING QA

Easily attain results for TG-142 recommended imaging QA, including:

- MV Imaging
- kV Imaging
- Cone-beam CT
- IGRT

### ● QUANTIFIABLE MACHINE QA

PIPSpro quickly analyzes specialized phantoms to provide precise, quantitative results for TG-142 machine QA procedures, including:

- MLC positioning (Static, IMRT), leakage and more
- Stereotactic (Winston-Lutz)
- Radiation Field / Light Field
- Star Shot analysis for rotational isocenter verification

### ● VALIDATED PERFORMANCE

Over 55 publications validate PIPSPRO and its specialized phantoms' performance, highlighting the software's compatibility with a variety of different imagers.

### ● REFINED USER-INTERFACE

New, streamlined user-interface simplifies QA testing. Easy-to-use tabs automate analysis, data storage, historical trending and reporting.

### ● CENTRALIZED DATABASE FOR TRENDING AND REPORTING

A central data hub allows for quick access to QA results from multiple computers, departments or sites.

# TG-142 Compliant QA Modules

## QC Module

Analysis of QC-3 and QcKV-1 Phantom images provides data for monthly MV/kV imaging tests, including:

- Spatial resolution\*
- Contrast-to-noise\*
- Overall noise\*

\*Meets TG-142 specification for comparison to baseline

Analysis of FC-2 Phantom images provides data for monthly mechanical and imaging QA procedures, including:

- Jaw position indicators (symmetric) – Meets TG-142 specification of measuring tolerances of 2mm
- Light/radiation field coincidence – Meets TG-142 specification of measuring tolerances of 2mm or 1% on a side
- Scaling, MV Imaging – Meets TG-142 specifications of measuring tolerances of  $\leq 2$ mm (non-SRS/SBRT) and  $\leq 1$ mm (SRS/SBRT)
- Scaling, kV imaging – Meets TG-142 specifications of measuring tolerances of  $\leq 2$ mm (non-SRS/SBRT) and  $\leq 1$ mm (SRS/SBRT)

## Star Shot Module

Analysis of star shot images provides data for annual mechanical procedures, including:

- Collimator rotation isocenter\*
- Gantry rotation isocenter\*
- Couch rotation isocenter\*

\*Meets TG-142 specification of measuring tolerances of  $\pm 1$ mm from baseline

## Stereotactic Module [COMPREHENSIVE]

Accurate, reproducible results:

- Automatically analyze EPID images of Winston-Lutz ball marker
- Assure patient position for the most precise stereotactic treatments using the calculated 3D offset

## MLC QA Module [COMPREHENSIVE]

Analysis of MLC Phantom images provides data for TG-142 procedures\*, including:

- Qualitative test [i.e., matched segments, aka 'picket fence']
- Setting vs. radiation field for two patterns [non-IMRT]
- Leaf position accuracy [IMRT]
- MLC Transmission [average of leaf and interleaf transmission, all energies]
- Leaf position repeatability
- Segmental IMRT [step and shoot] test

\*Exceeds TG-142 specifications with quantitative analysis

## IGRT Module [COMPREHENSIVE]

Track and trend results to prevent potential clinical issues:

- Imaging and treatment coordinate coincidence [kV, MV, CBCT] (single gantry angles and four cardinal angles)
- Positioning/repositioning [kV, MV, CBCT]

## CBCT Module [COMPREHENSIVE]

Automatic analysis of Catphan® phantoms provides data for monthly Cone-Beam CT QA Tests, including:

- Geometric distortion – Meets TG-142 specification of measuring tolerances of  $\leq 2$ mm (non-SRS/SBRT) and  $\leq 1$ mm (SRS/SBRT)
- Spatial Resolution\*
- Contrast\*
- HU Constancy\*
- Noise\*

\*Meets TG-142 specification for comparison to baseline

*The TG-58 Report encourages users to demand the use of a tool such as PIPspro QC Software at acceptance to help ensure that the EPID is indeed operating at or above specifications.\**

\* Michael G. Herman, et al, "Clinical Use of Electronic Portal Imaging: Report of AAPM Radiation Therapy Committee Task Group 58," **Med. Phys.** 28 (5), May 2001

The PIPspro QC Software (REF 91310) consists of the PIPspro QC Software, including modules for quality control, registration, image enhancement, and image analysis, plus the QC-3 and FC-2 Phantoms. The PIPspro Comprehensive Software (REF 91320) includes the MLC, QC-3, FC-2, and QcKV-1 Phantoms and additional modules for performing stereotactic radiosurgery QA, MLC QA, CBCT QA, and dewarping images.

**Over 55 publications validate the performance of PIPspro.**

## PIPspro Software SPECIFICATIONS

### OPERATING SYSTEM

MICROSOFT® WINDOWS® XP  
WINDOWS VISTA®  
WINDOWS 7

**PROCESSOR** Intel® or AMD®, 350 MHz or greater

**MEMORY** 64 MB (256 MB recommended)

**HARD DRIVE** 50 MB or greater

**SCREEN RESOLUTION** 800 x 600 (1024 x 768 recommended)

**CD-ROM DRIVE** 2X speed or greater

**SCREEN COLOR DEPTH** 256-bit or greater

**PRODUCT STANDARDS** Designed to meet IEC 60601-1-4 

**RECOMMENDED SOFTWARE** Microsoft Excel

Catphan® is a registered trademark of The Phantom Laboratory.  
Windows® is a registered trademark of Microsoft Corporation.  
Specifications subject to change without notice.

ORDERING  
INFORMATION

PIPSPRO SOFTWARE REF 91310